BioCUE BMA

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1. Hernigou, P., Poignard, A., Beaujean, F., Rouard, H. Percutaneous autologous bone-marrow grafting for nonunions. Influence of the number and concentration of progenitor cells. *Journal of Bone and Joint Surgery (American)*. 87(7): 1430–7, 2005.
2. Welch, Z.R., McKale, J.M., Woodell-May, J. Citrate-based Anticoagulant Does Not Hinder Stem Cell Viability and Concentration from Bone Marrow Aspirate. Submitted at 54th Annual Meeting of the Orthopaedic Research Society. March 2–5, 2008.
3. Hernigou, P., Mathieu, G., Poignard, A., Manicom, O., Beaujean, F., Rouard, H. Percutaneous autologous bone-marrow grafting for nonunions. Surgical technique. *Journal of Bone and Joint Surgery (American)*. 88 (Suppl 1 Pt 2): 322–7, 2006.
4. Siwach, R.C., Sangwan, S.S., Singh, R., Goel, A. Role of percutaneous bone marrow grafting in delayed unions, non-unions and poor regenerates. *Indian Journal of Medical Sciences*. 55(6): 326–36, 2001.
5. Aarvold, A., Smith, J.O., Tayton, E.R., Jones, A.M.H., Dawson, J.I., Lanham, S., Briscoe, A., Dunlop, D.G., and Oreffo, R.O.C. "A Tissue Engineering Strategy for the Treatment of Avascular Necrosis of the Femoral Head." *The Surgeon* 11.6 (2013): 319–25.
6. Martin, John R., Matthew T. Houdek, and Rafael J. Sierra. "Use of Concentrated Bone Marrow Aspirate and Platelet Rich Plasma during Minimally Invasive Decompression of the Femoral Head in the Treatment of Osteonecrosis." *Croatian Medical Journal* 54.3 (2013): 219–24.
7. Hernigou, P., Poignard, A., Zilber, S. and Rouard, H. "Cell Therapy of Hip Osteonecrosis with Autologous Bone Marrow Grafting." *Indian Journal of Orthopaedics* 43.1 (2009): 40.
8. Muschler, G.F., Matsukura, Y., Nitto, H., Boehm, C.A., Valdevit, A.D., Kambic, H.E., Davros, W.J., Easley, K.A., Powell, K.A. Selective retention of bone marrow-derived cells to enhance spinal fusion. *Clinical Orthopaedics and Related Research*. 432: 242–51, 2005.
9. Muschler, G.F., Nitto, H., Matsukura, Y., Boehm, C., Valdevit, A., Kambic, H., Davros, W., Powell, K., Easley, K. Spine fusion using cell matrix composites enriched in bone marrow-derived cells. *Clinical Orthopaedics and Related Research*. 407: 102–18, 2003.
10. Steinwachs, M.R., Waibl, B., Wopperer, S., Mumme, M. Matrix-associated chondroplasty: a novel platelet-rich plasma and concentrated nucleated bone marrow cell-enhanced cartilage restoration technique. *Arthroscopy Techniques*. 3(2) April: E279–e82, 2014.
11. Enea, D., Cecconi, S., Calcagno, S., Busilacchi, A., Manzotti, S., Kaps, C. and Gigante, A. "Single-stage Cartilage Repair in the Knee with Microfracture Covered with a Resorbable Polymer-based Matrix and Autologous Bone Marrow Concentrate." *The Knee* 20.6 (2013): 562–69.
12. Adachi, N., Ochi, M., Deie, M., Ito, Y. Transplant of mesenchymal stem cells and hydroxyapatite ceramics to treat severe osteochondral damage after septic arthritis of the knee. *Journal of Rheumatology*. 32(8):1615–18, 2005.
13. Wakitani, S., Imoto, K, Yamamoto, T., Saito, M., Murata, N., Yoneda, M. Human autologous culture expanded bone marrow mesenchymal cell transplantation for repair of cartilage defects in osteoarthritic knees. *Osteoarthritis Cartilage*. 10(3): 199–206, 2002.
14. Johnstone B, Yoo JU. Autologous mesenchymal progenitor cells in articular cartilage repair. *Clinical Orthopaedics and Related Research*. 367 (Suppl): S156–S162, 1999.
15. Jia, X., Peters, P. G. and Schon, L. "The Use of Platelet-Rich Plasma in the Management of Foot and Ankle Conditions." *Operative Techniques in Sports Medicine* 19.3 (2011): 177–84.
16. Adams, S. B., Lewis, J. S., Gupta, A. K., Parekh, S. G., Miller, S. D. and Schon, L. C. "Cannulated Screw Delivery of Bone Marrow Aspirate Concentrate to a Stress Fracture Nonunion: Technique Tip." *Foot & Ankle International* 34.5 (2013): 740–44.
17. Torres, J., Gutierrez, M., Lopes, M. A., Santos, J. D., Cabral, A. T., Pinto, R. and Van Eck, C. "Bone Marrow Stem Cells Added to a Hydroxyapatite Scaffold Result in Better Outcomes after Surgical Treatment of Intertrochanteric Hip Fractures." *BioMed Research International* (2014): 1–7. Bone Marrow Stem Cells Added to a Hydroxyapatite Scaffold Result in Better Outcomes after Surgical Treatment of Intertrochanteric Hip Fractures. Hindawi Publishing Corporation. Web. 06 Nov. 2014.
18. Tiedeman, J.J., Garvin, K.L., Kile, T.A., Connolly, J.F. The role of a composite, demineralized bone matrix and bone marrow in the treatment of osseous defects. *Orthopedics*. 18(12):1153–58, 1995.
19. Yamada, Y., Ueda, M., Naiki, T., Takahashi, M., Hata, K., Nagasaka, T. Autogenous injectable bone for regeneration with mesenchymal stem cells and platelet-rich plasma: tissue-engineered bone regeneration. *Tissue Engineering*. 10(5–6): 955–64, 2004.
20. Lucarelli, E., Fini, M., Beccheroni, A., Giavaresi, G., Di, B.C., Aldini, N.N. et al. Stromal stem cells and platelet-rich plasma improve bone allograft integration. *Clinical Orthopaedics and Related Research*. 435: 62–8, 2005.
21. Yamada, Y., Boo, J.S., Ozawa, R., Nagasaka, T., Okazaki, Y., Hata, K. et al. Bone regeneration following injection of mesenchymal stem cells and fibrin glue with a biodegradable scaffold. *Journal of Craniomaxillofacial Surgery*. 31(1): 27–33, 2003.
22. Yamada, Y., Ueda, M., Hibi, H., Baba, S. A novel approach to periodontal tissue regeneration with mesenchymal stem cells and platelet-rich plasma using tissue engineering technology: a clinical case report. *International Journal of Periodontics and Restorative Dentistry*. 26(4): 363–9, 2006.
23. Oyama, T., Nishimoto, S., Takeda, M. Alveolar bone regeneration utilizing b-TCP and platelet-rich plasma (PRP) derived from bone marrow aspirate. *Annals of Plastic Surgery*. 54(2): 222–3, 2005.
24. Yamada, Y., Ueda, M., Naiki, T., Nagasaka, T. Tissue-engineered injectable bone regeneration for osseointegrated dental implants. *Clinical Orthopaedics and Related Research*. 15(5): 589–97, 2004.

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Bibliografía (Cont.)

25. Docquier, P.L., Delloye, C. Treatment of aneurysmal bone cysts by introduction of demineralized bone and autogenous bone marrow. *Journal of Bone and Joint Surgery (American)*. 87(10): 2253–8, 2005.
26. Rougraff, B.T., Kling, T.J. Treatment of active unicameral bone cysts with percutaneous injection of demineralized bone matrix and autogenous bone marrow. *Journal of Bone and Joint Surgery (American)*. 84–A(6): 921–9, 2002.
27. Haynesworth, S.E., Kadiyala, S., Liang, L., Bruder, S.P. Mitogenic Stimulation of human mesenchymal stem cells by platelet releasate suggest a mechanism for enhancement of bone repair by platelet concentrates. *Transactions of the 48th Annual Meeting of the Orthopaedic Research Society*. 42: 0462, 2002.
28. Lucarelli, E., Beccheroni, A., Donati, D., Sangiorgi, L., Cenacchi, A., Del Vento A.M. et al. Platelet-derived growth factors enhance proliferation of human stromal stem cells. *Biomaterials*. 24(18): 3095–100, 2003.
29. Hakim, M., Grassmann, J. P., Betsch, M., Schneppendahl, J., Gehrman, S., Hakimi, A. R., Kropil, P., Sager, M., Herten, M., Wild, M., Windolf, J. and Jungbluth, P. "The Composite of Bone Marrow Concentrate and PRP as an Alternative to Autologous Bone Grafting." *PLOS ONE* E100143 9.6 (2014): 1-18. *PLOS ONE*. June 2014. Web. 06 Nov. 2014.

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Biomet, Inc.
Biomet UK, Ltd.